

OPHTHALMIC OPERATION TABLE

MO-2012 SS



Two M Ophthotronics

Hi-Tech Ophtho Equipment

OPHTHALMIC OPERATION TABLE MO-2012



| TECHNICAL SPECIFICATIONS | |
|--------------------------|------------------------------------|
| Height range min./ max | 640 mm / 840 mm |
| Lifting speed | 8 mm/s |
| Total length | 1880 mm (with head) |
| Total width | 600 mm |
| Base size | 920 mm / 540 mm |
| Input Power Supply | 230 V 50/60 Hz. 110 V 50/60 Hz. |
| Weight capacity | 170 Kgs. |
| Gross weight | 130 Kgs. |

[Please Note: Specification are subject to modifications, as per emerging market demands.]

Motorised Vertical Height Adjustment

enables the ophthalmologist to perform from a convenient, comfortable standing or sitting position.

Optimally Positioned Head Rest

with manually controlled up & down motion option allows the surgeon easy access to the patient, while offering the patient comfortable head positioning. Plus, there is ample leg room, and sufficient space for equipment foot pedals.

Noiseless Dc Motor-driven

Adjustable Wrist Support

Offers surgeons convenient stability during operations.

Built-in Head Low and Head Up Position ensures safety during emergencies.



Two M Ophthotronics
Hi-Tech Ophtho Equipment

www.twom.co.in

Regd. Off.: No.-35, Satnam Industrial Estate, Nagarvell Hanuman Road,
Near SBI Bank, Amraiwadi, Ahmedabad-380026. Gujarat, INDIA.
E-mail : twomophtho@gmail.com

Branch Off.: Unit No. 110, Sahar Cargo Estate, V.M. Shah Marg.
Nr. Tarun Bharat Society, Chakala, Andheri (E), Mumbai-400099, INDIA.
Ph.: 022 2832 5213 | **Fax** : 022 2825 0347
E-mail : twom.raj@gmail.com