

OPHTHALMIC WORKSTATION

MR-2012



**TWO-M OPHTHOTRONICS
PVT. LTD.**

www.twom.in

"We are Digital"

CHAIR UNIT

Fully comfortable rotating chair,
 Silent CE marked DC motor
 gives smooth up-down & back-forward
 movement, Adjustable footrest
 which can be fully raised to
 sleeping position, Foot switch for
 chair up/down movement control.
**(Available in Right hand (Standard)
 & Left Hand models)**




Trail Set Storage



Writing Table With Mini Storage

OPHTHALMIC WORKSTATION MR-2012



TECHNICAL SPECIFICATIONS

Chair height Min. /Max.	560 / 760mm
Chair Elevation Stroke	200 mm
Instrument Table Height	875 mm
Maximum Area Required	1450 x 1250 mm in erect position (Without Writing Pad)
Back-rest	Motorized
Table Top Size	32" x 16"
Rotation	180 circular
Load Lifting Capacity	175 kgs.
Auxiliary Output	230v 50 Hz/ 4 connection
Voltage on table	3v/6v/7.5v/9v/12v/15v
Power Supply	230 v 50/60 Hz 5 Amp.
Gross Weight	240 Kgs.

Please Note: Specification are subject to modifications,
 as per emerging market demands.



**TWO-M OPHTHOTRONICS
 PVT. LTD.**

Hi-Tech Ophtho Equipment

www.twom.in

Regd. Off.: No.-35, Satnam Industrial Estate, Nagarvell Hanuman Road,
 Near SBI Bank, Amraiwadi, Ahmedabad-380026. Gujarat, INDIA.

E-mail : twomophtho@gmail.com

Branch Off.: Unit No. 110, Sahar Cargo Estate, V.M. Shah Marg.
 Nr. Tarun Bharat Society, Chakala, Andheri (E), Mumbai-400099, INDIA.

Ph.: 022 2832 5213 | **E-mail :** twom.raj@gmail.com