



**TWO-M OPHTHOTRONICS  
PVT. LTD.**



**GENERAL**  
EXAMINATION CHAIR  
ME-2012

Fully Flexible

Rotating Recliner Chair

High Stability



**TWO-M OPHTHOTRONICS  
PVT. LTD.**

*"We are Digital"*



## TWO-M OPHTHOTRONICS PVT. LTD.

### Tehchnical Specification

- Fully Flexible, Rotating Recliner Chair
- Noiseless CE-Marked, Imported Dc Motor-driven
- Up & Down Movements And Inclination
- High Stability
- Comfortable Seat.
- Arm Rest With Swing-up Option
- Adjustable Footrest, Which Can Be Fully Raised
- To Sleeping Position
- Ideal Movement Control From Chair Side Top
- Foot Switch For Up / Down Movement Control



### Specification

### ME-2012

Height min. / max.	560mm / 760mm $\pm$ 10mm
Maximum load	175kgs.
Backrest	Electrical & Manual
Rotation	180-degree
Power	300 VA
Input Power Supply	230 V 50/60 Hz. 110 V 50/60 Hz.
Total weight	100 Kgs.



## TWO-M OPHTHOTRONICS PVT. LTD.

Hi-Tech Ophtho Equipment

[www.twom.in](http://www.twom.in)

Regd. Off.: No.-35, Satnam Industrial Estate, Nagarvell Hanuman Road,  
Near SBI Bank, Amraiwadi, Ahmedabad-380026. Gujarat, INDIA.

**E-mail :** [twomophtho@gmail.com](mailto:twomophtho@gmail.com)

Branch Off.: Unit No. 110, Sahar Cargo Estate, V.M. Shah Marg.  
Nr. Tarun Bharat Society, Chakala, Andheri (E), Mumbai-400099, INDIA.

**Ph.:** 022 2832 5213 | **E-mail :** [twom.raj@gmail.com](mailto:twom.raj@gmail.com)