



**TWO-M OPHTHOTRONICS
PVT. LTD.**



Model - DV003 Intergrated Camera Module **Third-I**

FEAUTRES

- Image Capturing Device with Single-CCD Camera Module & Control Unit Built - in, For Exclusive Use with Operating Microscopes & slitlamps.
- No separate beam splitter or camera adapter is rquires.
- The power source & singal are transmitted on single cable.
- User Name on Display
- Special Port to Adjust the Colour, Brightness etc
- Shutter illumination can be Adjusted

TECHNICAL SPECIFIACATIONS

CHIP	SONY (SINGLE CHIP)
SENSOR SIZE	1/3" CCD SENSOR
RESOLUTION	2 MP
WHITE BALANCE	AUTO WHITE TRACE
GAIN	AUTO GAIN CONTROL
TV SYSTEM	PAL
POWER SUPPLY	12V DC 500MA
DIMENSION	91 X 70 X 52 (MM)
WEIGHT	410 GRAMS



**TWO-M OPHTHOTRONICS
PVT. LTD.**

Hi-Tech Ophtho Equipment

www.twom.in

Regd. Off.: No.-35, Satnam Industrial Estate, Nagarvell Hanuman Road,
Near SBI Bank, Amraiwadi, Ahmedabad-380026. Gujarat, INDIA.

E-mail : twomophtho@gmail.com

Branch Off.: Unit No. 110, Sahar Cargo Estate, V.M. Shah Marg.
Nr. Tarun Bharat Society, Chakala, Andheri (E), Mumbai-400099, INDIA.

Ph.: 022 2832 5213 | **E-mail :** twom.raj@gmail.com