

# *DIGITAL* REFRACTION UNIT

## **KEPLER-17**

*LED TOUCH PANEL*


 *BLUE TOOTH*

*AUTOMATED ARM*



**TWO-M OPHTHOTRONICS  
PVT. LTD.**

*“We are Digital”*



**Automated Arm**

# DIGITAL REFRACTION UNIT

## KEPLER-17



**LED touch panel**




**Luxurious leather finish**



**Multifunction footswitch  
(controls patient chair & phoropter arm)**



**control the unit through "APP"**



### TECHNICAL SPECIFICATIONS

Chair height Min./Max.	560 / 760mm
Chair Elevation Stroke	200 mm
Instrument Table Height	875 mm
Maximum Area Required	2120 x 1700 mm in erect position
Back-rest	Motorized
Table Top Size	32" x 16"
Rotation	180 circular
Load Lifting Capacity	175 kgs.
Auxiliary Output	230v 50 Hz/ 4 connection
Voltage on table	6v/12v/15v
Power Supply	230 v / 110 v 50/60 Hz
Gross Weight	260 Kgs.

Please Note: Specification are subject to modifications, as per emerging market demands.



**TWO-M OPHTHOTRONICS  
PVT. LTD.**

Hi-Tech Ophtho Equipment

[www.twom.in](http://www.twom.in)

Regd. Off.: No.-35, Satnam Industrial Estate, Nagarvell Hanuman Road,  
Near SBI Bank, Amraiwadi, Ahmedabad-380026. Gujarat, INDIA.

**E-mail :** [twomophtho@gmail.com](mailto:twomophtho@gmail.com)

Branch Off.: Unit No. 110, Sahar Cargo Estate, V.M. Shah Marg.  
Nr. Tarun Bharat Society, Chakala, Andheri (E), Mumbai-400099, INDIA.

**Ph.:** 022 2832 5213 | **E-mail :** [twom.raj@gmail.com](mailto:twom.raj@gmail.com)