

OPHTHALMIC WORKSTATION

MR-2012 OPTICAL



**TWO-M OPHTHOTRONICS
PVT. LTD.**

"We are Digital"

www.twom.in

Optional Instrument

Visual Acuity Chart
Model. EYE-Q




REMOTE CONTROLLED
VISION DRUM

Model : MVD2015



OPHTHALMIC WORKSTATION MR-2012 OPTICAL



TECHNICAL SPECIFICATIONS

Chair height Min. /Max.	560 / 760mm
Chair Elevation Stroke	200 mm
Instrument Table Height	875 mm
Maximum Area Required	1450 x 1250 mm in erect position
Back-rest	Manual
Table Top Size	32" x 16"
Rotation	180 circular
Load Lifting Capacity	175 kgs.
Auxiliary Output	230v 50 Hz/ 4 connection
Voltage on table	3v/6v/7.5v/9v/12v/15v
Power Supply	230 v 50/60 Hz 5 Amp.
Gross Weight	235 Kgs.

Please Note: Specification are subject to modifications,
as per emerging market demands.



**TWO-M OPHTHOTRONICS
PVT. LTD.**

Hi-Tech Ophtho Equipment

www.twom.in

Regd. Off.: No.-35, Satnam Industrial Estate, Nagarvell Hanuman Road,
Near SBI Bank, Amraiwadi, Ahmedabad-380026. Gujarat, INDIA.

E-mail : twomophtho@gmail.com

Branch Off.: Unit No. 110, Sahar Cargo Estate, V.M. Shah Marg.
Nr. Tarun Bharat Society, Chakala, Andheri (E), Mumbai-400099, INDIA.

Ph.: 022 2832 5213 | **E-mail :** twom.raj@gmail.com