

Two M Ophthotronics

Hi-tech Ophtho Equipment

AN ISO 9001 : 2008 CERTIFIED COMPANY



Pot-360^{Gold}
OPHTHALMIC
MOBILE OPERATING TABLE
STAINLESS STEEL

POT-360Gold

OPHTHALMIC

MOBILE OPERATING TABLE

STAINLESS STEEL



- POT-360Gold Stainless steel mobile operating table is meets all the requirements of ophthalmologists.
- The mobile operating table POT-360Gold offers many extra features and a large range of positioning possibilities.
- The special design of the castors enables both straight-line travel as well as free-wheeling

FEAUTRES

- EASY TRANSPORTATION OF PATIENTS
- TRENDELENBURG / REVERSE TRENDELENBURG
- MORE LEG SPACE
- IMPORTED DC MOTOR DRIVEN
- FOOT SWITCH & HAND CONTROL OPERATED
- BATTERY BACK-UP(OPTIONAL)

TECHNICAL SPECIFIACATIONS

HEIGHT RANGE MIN. / MAX	680MM / 880MM
LIFTING SPEED	8MM / S
TOTAL LENGTH	1880MM (WITH HEAD)
TOTAL WIDTH	600MM
BASE SIZE	840MM / 520MM
INPUT POWER SUPPLY	230 V 50/60 HZ.
	110 V 50/60 HZ.
WEIGHT LIFTING CAPACITY	250 KGS.
GROSS WEIGHT	130KGS.

NIRMTCOS 01/15 : 99098 95777 | 97276 57002



Two M Ophthotronics

Hi-tech Ophtho Equipment
www.twom.co.in

Regd. off. : No-35,satnam Industrial Estate,
Nagarvell Hanuman Road, Near SBI Bank, Amraiwadi,
Ahmedabad-380026. Gujarat, India.
E-mail : twomophtho@gmail.com

Branch Off. : Unit No. 110, Sahar Cargo Estate, V.M. Shah Marg,
Near Tarun Bharat Society, Chakala, Andheri (E),
Mumbai - 400 099, India.
Ph.: +91 22 2832 5213, Fax : +91 22 2825 0347
E-mail : twom.raj@gmail.com