



**TWO-M  
OPHTHOTRONICS  
PVT. LTD.**



**MANUAL  
PHOROPTER  
ARM**

**PAM-25**

MANUAL  
PHOROPTER  
ARM

**PAM-25**

Specification

Main Holding Collar  
**42mm**


Refractor Adjustment  
**70mm**

Up/Down (Elevation)  
**280mm**

Refractor Mounting Pin  
**15.7mm**

Load capacity  
**5kg**

Gross Weight  
**10kg**



**TWO-M OPHTHOTRONICS  
PVT. LTD.**  
Hi-Tech Ophtho Equipment



CORPORATE OFFICE

101/102, Sahar Cargo Estate, V.M Shah Marg,  
Chakala, Andheri (East), Mumbai - 400 099, India.  
Mo.: +91 9820351012 | +91 85913 77382  
E : twom.mumbai@gmail.com

FACTORY ADDRESS

35, Satnam Industrial Estate, Nagarvel Hanuman  
Road, Amraiwadi, Ahmedabad - 380 026, India.  
Mo.: +91 99786 02571 | +91 87782 17424  
E : twomamd@gmail.com | twomaccamd@gmail.com