



**TWO-M
OPHTHOTRONICS
PVT. LTD.**

OPHTHALMIC WORKSTATION

MR-2012
VINTAGE



OPHTHALMIC WORKSTATION

MR-2012 VINTAGE



OPEN TRAY FOR
TRIAL CASE



Technical Specifications

Chair Height Min./max.
560/760mm

Chair Elevation Stroke
200mm

Load Lifting Capacity
175 Kg

Rotation
180 Circular

Power Supply
230v 50/60hz 5amp

Gross Weight
200kg


Instrument Table Height
875mm

Maximum Area Required
1970*1420mm

Table To Size
860*400mm

Tray Size
530*330 mm

Recline
Manual



**TWO-M OPHTHOTRONICS
PVT. LTD.**
Hi-Tech Ophtho Equipment

CORPORATE OFFICE

101/102, Sahar Cargo Estate, V.M Shah Marg,
Chakala, Andheri (East), Mumbai - 400 099, India.
Mo.: +91 9820351012 | +91 85913 77382
E : twom.mumbai@gmail.com

FACTORY ADDRESS

35, Satnam Industrial Estate, Nagarvel Hanuman
Road, Amraiwadi, Ahmedabad - 380 026, India.
Mo.: +91 99786 02571 | +91 87782 17424
E : twomamd@gmail.com | twomaccamd@gmail.com

